

3M Skin & Wound Care Division
Medical Specialties



Trusted.
Innovative.
Global.

Customer Focused Solutions
for Medical Business to Business Markets

New Silicone Pressure Sensitive Adhesive Tapes

Introducing a new skin-friendly, adhesive technology for use in the design and manufacturing of medical devices. With new, proprietary technology, we can apply silicone adhesives to multiple substrates to achieve a variety of properties and configurations. This new line of silicone pressure sensitive tapes will be manufactured with the consistency, quality and scale that 3M has exhibited over its 60 plus years in the Healthcare adhesives market.



3M™ Silicone Adhesive Tape Portfolio

Features:

- 4 mils silicone-based pressure sensitive adhesive
- Standard or premium release film liners available
- EtO sterilization compatible
- Wide range of substrates/backings available

Typical Benefits of Using Silicone Adhesives*:

- Repositionable
- Gentle adhesion to skin
- Minimizes skin and hair trauma during removal
- Conformable
- Hypoallergenic
- Consistent adhesion to skin during wear
- Latex-free

Why Select 3M Silicone Adhesive Tapes?*

- Excellent tack and adhesion properties
- Strong bond between adhesive and tape backing
- Ability to coat adhesive on a range of tape backings
- Tin, platinum and peroxide-free silicone technology
- 3M quality assurance and high level of support
- Clinical summaries

Catalog Number	Description
2475	3M™ 0.6-mil Thermoplastic Elastomer (TPE) Tape with Silicone Adhesive
2475P	3M™ 0.6-mil TPE Tape with Silicone Adhesive on Premium Release Liner
2476	3M™ Spunlace Polyester Nonwoven Tape with Silicone Adhesive
2476P	3M™ Spunlace Polyester Nonwoven Tape with Silicone Adhesive on Premium Release Liner
2477	3M™ Double-Coated TPE Tape with Silicone and Acrylate Adhesives
2477P	3M™ Double-Coated TPE Tape with Silicone and Acrylate Adhesives on Premium Release Liner

Ordering Information

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: www.3M.com/medicalspecialties for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative.

Product and Safety Information: User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

Warranty Information: All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of Remedies: If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, indirect, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict liability.

* Prototype typical properties. Properties listed based on very limited data and subject to change.



Medical Specialties

3M Center, Building 275-5W-05

St. Paul, MN 55144-1000

USA

1-800-228-3957

www.3M.com/medicalspecialties

3M is a trademark of 3M.

Please recycle. Printed in U.S.A.

© 2012, 3M. All rights reserved.

70-2010-9010-0

